

Dextrin Conjugated Antibody

Catalog No: #C56467

Package Size: #C56467-AF350 100ul #C56467-AF405 100ul #C56467-AF488 100ul #C56467-AF555 100ul #C56467-AF594 100ul
#C56467-AF647 100ul #C56467-AF680 100ul #C56467-AF750 100ul #C56467-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dextrin Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Dextrin Antibody detects endogenous levels of total Dextrin
Immunogen Description	A synthesized peptide derived from human Dextrin
Conjugates	AF Dye
Other Names	ACTDP; ADF; corn1; Dextrin; DSN; Dstn; HEL32; Sid 23; sid23p;
Accession No.	Uniprot:P60981
Calculated MW	17kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.