

## HSC70 Interacting Protein HIP Conjugated Antibody

Catalog No: #C56475

Package Size: #C56475-AF350 100ul #C56475-AF405 100ul #C56475-AF488 100ul #C56475-AF555 100ul #C56475-AF594 100ul  
#C56475-AF647 100ul #C56475-AF680 100ul #C56475-AF750 100ul #C56475-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                         |
|-----------------------|---------------------------------------------------------------------------------------------------------|
| Product Name          | HSC70 Interacting Protein HIP Conjugated Antibody                                                       |
| Host Species          | Rabbit                                                                                                  |
| Clonality             | Monoclonal                                                                                              |
| Isotype               | Rabbit IgG                                                                                              |
| Purification          | Affinity-chromatography                                                                                 |
| Species Reactivity    | Human Mouse                                                                                             |
| Specificity           | HSC70 Interacting Protein HIP Antibody detects endogenous levels of total HSC70 Interacting Protein HIP |
| Immunogen Description | A synthesized peptide derived from human HSC70 Interacting Protein HIP                                  |
| Conjugates            | AF Dye                                                                                                  |
| Other Names           | HIP; HOP; P48; AAG2; SNC6; HSPABP; FAM10A1; FAM10A4; HSPABP1; ST13;                                     |
| Accession No.         | Uniprot:P50502                                                                                          |
| Calculated MW         | 47kDa                                                                                                   |
| Storage               | Store at 4°C in dark for 6 months                                                                       |

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.