

# Trefoil Factor 3 Conjugated Antibody

Catalog No: #C56482

Package Size: #C56482-AF350 100ul #C56482-AF405 100ul #C56482-AF488 100ul #C56482-AF555 100ul #C56482-AF594 100ul  
#C56482-AF647 100ul #C56482-AF680 100ul #C56482-AF750 100ul #C56482-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                               |
|-----------------------|-------------------------------------------------------------------------------|
| Product Name          | Trefoil Factor 3 Conjugated Antibody                                          |
| Host Species          | Rabbit                                                                        |
| Clonality             | Monoclonal                                                                    |
| Isotype               | Rabbit IgG                                                                    |
| Purification          | Affinity-chromatography                                                       |
| Species Reactivity    | Human                                                                         |
| Specificity           | Trefoil Factor 3 Antibody detects endogenous levels of total Trefoil Factor 3 |
| Immunogen Description | A synthesized peptide derived from human Trefoil Factor 3                     |
| Conjugates            | AF Dye                                                                        |
| Other Names           | hITF; ITF; mITF; P1B; TFF3; TFI; TREFOIL;                                     |
| Accession No.         | Uniprot:Q07654                                                                |
| Calculated MW         | 10kDa                                                                         |
| Storage               | Store at 4°C in dark for 6 months                                             |

## Application Details

### Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.