

MARK3 Conjugated Antibody

Catalog No: #C56483

Package Size: #C56483-AF350 100ul #C56483-AF405 100ul #C56483-AF488 100ul #C56483-AF555 100ul #C56483-AF594 100ul
#C56483-AF647 100ul #C56483-AF680 100ul #C56483-AF750 100ul #C56483-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MARK3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Mouse, Human
Specificity	MARK3 Antibody detects endogenous levels of total MARK3
Immunogen Description	A synthesized peptide derived from human MARK3
Conjugates	AF Dye
Other Names	C-TAK1; cTAK1; Emk2; KP78; Mark3; PAR1A;
Accession No.	Uniprot:P27448
Calculated MW	86kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.