

Vav1/2/3 Conjugated Antibody

Catalog No: #C56485

Package Size: #C56485-AF350 100ul #C56485-AF405 100ul #C56485-AF488 100ul #C56485-AF555 100ul #C56485-AF594 100ul
#C56485-AF647 100ul #C56485-AF680 100ul #C56485-AF750 100ul #C56485-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vav1/2/3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	Vav1/2/3 Antibody detects endogenous levels of total Vav1/2/3
Immunogen Description	A synthesized peptide derived from human Vav1/2/3
Conjugates	AF Dye
Other Names	p95; p95vav; Protein vav 1; Protein vav 2; Protein vav 3; VAV; vav-T;
Accession No.	Uniprot:P15498/P52735/Q9UKW4
Calculated MW	98kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.