

METTL3 Conjugated Antibody

Catalog No: #C56489

Package Size: #C56489-AF350 100ul #C56489-AF405 100ul #C56489-AF488 100ul #C56489-AF555 100ul #C56489-AF594 100ul
#C56489-AF647 100ul #C56489-AF680 100ul #C56489-AF750 100ul #C56489-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	METTL3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	METTL3 Antibody detects endogenous levels of total METTL3
Immunogen Description	A synthesized peptide derived from human METTL3
Conjugates	AF Dye
Other Names	IME4; M6A; Methyltransferase like protein 3; METTL3; MTA70; N6 adenosine methyltransferase 70 kDa subunit;
Accession No.	Uniprot:Q86U44
Calculated MW	70kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.