

# GluA4 Conjugated Antibody

Catalog No: #C56497



Package Size: #C56497-AF350 100ul #C56497-AF405 100ul #C56497-AF488 100ul #C56497-AF555 100ul #C56497-AF594 100ul  
#C56497-AF647 100ul #C56497-AF680 100ul #C56497-AF750 100ul #C56497-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                         |
|-----------------------|---------------------------------------------------------|
| Product Name          | GluA4 Conjugated Antibody                               |
| Host Species          | Rabbit                                                  |
| Clonality             | Monoclonal                                              |
| Isotype               | Rabbit IgG                                              |
| Purification          | Affinity-chromatography                                 |
| Species Reactivity    | Human Mouse Rat                                         |
| Specificity           | GluA4 Antibody detects endogenous levels of total GluA4 |
| Immunogen Description | A synthesized peptide derived from human GluA4          |
| Conjugates            | AF Dye                                                  |
| Other Names           | AMPA4; GluA4; GluR4; GLUR4C; GLURD; Gria4;              |
| Accession No.         | Uniprot:P48058                                          |
| Calculated MW         | 101kDa                                                  |
| Storage               | Store at 4°C in dark for 6 months                       |

## Application Details

Suggested Dilution:

- AF350 conjugated: most applications: 1: 50 - 1: 250
- AF405 conjugated: most applications: 1: 50 - 1: 250
- AF488 conjugated: most applications: 1: 50 - 1: 250
- AF555 conjugated: most applications: 1: 50 - 1: 250
- AF594 conjugated: most applications: 1: 50 - 1: 250
- AF647 conjugated: most applications: 1: 50 - 1: 250
- AF680 conjugated: most applications: 1: 50 - 1: 250
- AF750 conjugated: most applications: 1: 50 - 1: 250
- Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.