

Syntaxin Conjugated Antibody

Catalog No: #C56510

Package Size: #C56510-AF350 100ul #C56510-AF405 100ul #C56510-AF488 100ul #C56510-AF555 100ul #C56510-AF594 100ul
#C56510-AF647 100ul #C56510-AF680 100ul #C56510-AF750 100ul #C56510-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Syntaxin Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Syntaxin Antibody detects endogenous levels of total Syntaxin
Immunogen Description	A synthesized peptide derived from human Syntaxin
Conjugates	AF Dye
Other Names	EPIM; Epimorphin; EPM; HPC 1; P35 1; STX1A; Stx1b; STX1B1; STX1B2; STX2; STX2A; STX2B; STX2C; STX3; STX3A;
Accession No.	Uniprot:P61266
Calculated MW	33kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.