

AKR1C3 Conjugated Antibody

Catalog No: #C56513



Package Size: #C56513-AF350 100ul #C56513-AF405 100ul #C56513-AF488 100ul #C56513-AF555 100ul #C56513-AF594 100ul
#C56513-AF647 100ul #C56513-AF680 100ul #C56513-AF750 100ul #C56513-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	AKR1C3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	AKR1C3 Antibody detects endogenous levels of total AKR1C3
Immunogen Description	A synthesized peptide derived from human AKR1C3
Conjugates	AF Dye
Other Names	DD3; DDX; HA1753; HAKRB; HAKRe; HSD17B5; KIAA0119; hluPGFS;
Accession No.	Uniprot:P42330
Calculated MW	37kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB:1:500~1:1000
ICC/IF:1:50~1:200
IP:1:30
FC:1:50

Note: This product is for in vitro research use only and is not intended for use in humans or animals.