

COPS3 Conjugated Antibody

Catalog No: #C56528

Package Size: #C56528-AF350 100ul #C56528-AF405 100ul #C56528-AF488 100ul #C56528-AF555 100ul #C56528-AF594 100ul
#C56528-AF647 100ul #C56528-AF680 100ul #C56528-AF750 100ul #C56528-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	COPS3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	COPS3 Antibody detects endogenous levels of total COPS3
Immunogen Description	A synthesized peptide derived from human COPS3
Conjugates	AF Dye
Other Names	COP9 subunit 3; cops3; CSN3; JAB 1 containing signalosome subunit 3; SGN3; Signalosome subunit 3;
Accession No.	Uniprot:Q9UNS2
Calculated MW	48 kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.