

CD16 Conjugated Antibody

Catalog No: #C56554

Package Size: #C56554-AF350 100ul #C56554-AF405 100ul #C56554-AF488 100ul #C56554-AF555 100ul #C56554-AF594 100ul
#C56554-AF647 100ul #C56554-AF680 100ul #C56554-AF750 100ul #C56554-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD16 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR1704
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	CD16 Antibody detects endogenous levels of total CD16
Immunogen Description	A synthesized peptide derived from human CD16
Conjugates	AF Dye
Other Names	IGFR3; CD16; CD16a; IMD20;
Accession No.	Uniprot:P08637
Calculated MW	42kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.