

SQSTM1 / p62 Conjugated Antibody

Catalog No: #C56556

Package Size: #C56556-AF350 100ul #C56556-AF405 100ul #C56556-AF488 100ul #C56556-AF555 100ul #C56556-AF594 100ul
#C56556-AF647 100ul #C56556-AF680 100ul #C56556-AF750 100ul #C56556-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SQSTM1 / p62 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	SQSTM1 / p62 Antibody detects endogenous levels of total SQSTM1 / p62
Immunogen Description	A synthesized peptide derived from human SQSTM1 / p62
Conjugates	AF Dye
Other Names	A170; DMRV; EBIAP; ORCA; OSF-6; Osi; OSIL; p60; p62; p62B; PDB3; Sequestosome 1; Sqstm1; STAP; STONE14; Ubiquitin binding protein p62; ZIP; ZIP3;
Accession No.	Uniprot:Q13501
Calculated MW	62 kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.