

TrkB Conjugated Antibody

Catalog No: #C56557

Package Size: #C56557-AF350 100ul #C56557-AF405 100ul #C56557-AF488 100ul #C56557-AF555 100ul #C56557-AF594 100ul
#C56557-AF647 100ul #C56557-AF680 100ul #C56557-AF750 100ul #C56557-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TrkB Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Mouse Rat
Specificity	TrkB Antibody detects endogenous levels of total TrkB
Immunogen Description	A synthesized peptide derived from human TrkB
Conjugates	AF Dye
Other Names	BDNF tropomyosine receptor kinase B; GP145 TrkB; GP145-TrkB/GP95-TrkB; GP95 TrkB; Neurotrophic receptor tyrosine kinase 2; Ntrk2; RATTRKB1; Tkrb; TRKB; TrkB tyrosine kinase; TRKB1; Tropomyosin related kinase B;
Accession No.	Uniprot:P15209
Calculated MW	90,140kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.