

TCP1 beta Conjugated Antibody

Catalog No: #C56563

Package Size: #C56563-AF350 100ul #C56563-AF405 100ul #C56563-AF488 100ul #C56563-AF555 100ul #C56563-AF594 100ul
#C56563-AF647 100ul #C56563-AF680 100ul #C56563-AF750 100ul #C56563-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TCP1 beta Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	TCP1 beta Antibody detects endogenous levels of total TCP1 beta
Immunogen Description	A synthesized peptide derived from human TCP1 beta
Conjugates	AF Dye
Other Names	CCTB; 99D8.1; PRO1633; CCT-beta; MGC142074; MGC142076; TCP-1-beta;
Accession No.	Uniprot:P78371
Calculated MW	57kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.