

PSMA2 Conjugated Antibody

Catalog No: #C56577

Package Size: #C56577-AF350 100ul #C56577-AF405 100ul #C56577-AF488 100ul #C56577-AF555 100ul #C56577-AF594 100ul
#C56577-AF647 100ul #C56577-AF680 100ul #C56577-AF750 100ul #C56577-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PSMA2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	PSMA2 Antibody detects endogenous levels of total PSMA2
Immunogen Description	A synthesized peptide derived from human PSMA2
Conjugates	AF Dye
Other Names	HC3; PMSA2; Proteasome component C3; PSC2; PSC3; psmA2;
Accession No.	Uniprot:P25787
Calculated MW	25kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.