

BCL2L12 Conjugated Antibody

Catalog No: #C56580

Package Size: #C56580-AF350 100ul #C56580-AF405 100ul #C56580-AF488 100ul #C56580-AF555 100ul #C56580-AF594 100ul
#C56580-AF647 100ul #C56580-AF680 100ul #C56580-AF750 100ul #C56580-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BCL2L12 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	BCL2L12 Antibody detects endogenous levels of total BCL2L12
Immunogen Description	A synthesized peptide derived from human BCL2L12
Conjugates	AF Dye
Other Names	BCL2 like 12; BCL2L12; BPR;
Accession No.	Uniprot:Q9HB09
Calculated MW	32kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.