

Calpain S1 Conjugated Antibody

Catalog No: #C56606

Package Size: #C56606-AF350 100ul #C56606-AF405 100ul #C56606-AF488 100ul #C56606-AF555 100ul #C56606-AF594 100ul #C56606-AF647 100ul #C56606-AF680 100ul #C56606-AF750 100ul #C56606-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Calpain S1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Calpain S1 Antibody detects endogenous levels of total Calpain S1
Immunogen Description	A synthesized peptide derived from human Calpain S1
Conjugates	AF Dye
Other Names	Calpain4; CalpainS1; CANP; CANPS; CAPN 4; CAPN S1; CAPN4; CAPNS 1; CAPNS; Capns1; CDPS; CSS1;
Accession No.	Uniprot:P04632
Calculated MW	28kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.