

Cyclophilin 40 Conjugated Antibody

Catalog No: #C56619

Package Size: #C56619-AF350 100ul #C56619-AF405 100ul #C56619-AF488 100ul #C56619-AF555 100ul #C56619-AF594 100ul
#C56619-AF647 100ul #C56619-AF680 100ul #C56619-AF750 100ul #C56619-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclophilin 40 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Cyclophilin 40 Antibody detects endogenous levels of total Cyclophilin 40
Immunogen Description	A synthesized peptide derived from human Cyclophilin 40
Conjugates	AF Dye
Other Names	Cyclophilin40; CyclophilinD; CYP40; CYPD; PPlase; PPlase D; Ppid; Rotamase; Rotamase D;
Accession No.	Uniprot:Q08752
Calculated MW	41kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.