

RAP1A + RAP1B Conjugated Antibody

Catalog No: #C56620

Package Size: #C56620-AF350 100ul #C56620-AF405 100ul #C56620-AF488 100ul #C56620-AF555 100ul #C56620-AF594 100ul
#C56620-AF647 100ul #C56620-AF680 100ul #C56620-AF750 100ul #C56620-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RAP1A + RAP1B Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	RAP1A + RAP1B Antibody detects endogenous levels of total RAP1A + RAP1B
Immunogen Description	A synthesized peptide derived from human RAP1A + RAP1B
Conjugates	AF Dye
Other Names	C21KG; G 22K; K REV; KREV 1; KREV; KREV1; RAL1B; RAP 1; RAP 1A; RAP 1B; RAP1; RAP1B; SMGP21;
Accession No.	Uniprot:P61224/P62834
Calculated MW	21kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.