

CAD Conjugated Antibody

Catalog No: #C56656

Package Size: #C56656-AF350 100ul #C56656-AF405 100ul #C56656-AF488 100ul #C56656-AF555 100ul #C56656-AF594 100ul
#C56656-AF647 100ul #C56656-AF680 100ul #C56656-AF750 100ul #C56656-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CAD Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	CAD Antibody detects endogenous levels of total CAD
Immunogen Description	A synthesized peptide derived from human CAD
Conjugates	AF Dye
Other Names	Aspartate transcarbamylase; CAD protein; CAD trifunctional protein; Carbamoyl phosphate synthetase 2; CPSase ATCase DHOase; Dihydroorotase;
Accession No.	Uniprot:P27708
Calculated MW	260kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.