

# Granulin Conjugated Antibody

Catalog No: #C56657

Package Size: #C56657-AF350 100ul #C56657-AF405 100ul #C56657-AF488 100ul #C56657-AF555 100ul #C56657-AF594 100ul  
#C56657-AF647 100ul #C56657-AF680 100ul #C56657-AF750 100ul #C56657-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                             |
|-----------------------|---------------------------------------------------------------------------------------------|
| Product Name          | Granulin Conjugated Antibody                                                                |
| Host Species          | Rabbit                                                                                      |
| Clonality             | Monoclonal                                                                                  |
| Isotype               | Rabbit IgG                                                                                  |
| Purification          | Affinity-chromatography                                                                     |
| Species Reactivity    | Human                                                                                       |
| Specificity           | Granulin Antibody detects endogenous levels of total Granulin                               |
| Immunogen Description | A synthesized peptide derived from human Granulin                                           |
| Conjugates            | AF Dye                                                                                      |
| Other Names           | Acrogranin; CLN11; GEP; GP88; Granulins; GRN; PCDGF; PEPI; PGRN; Proepithelin; Progranulin; |
| Accession No.         | Uniprot:P28799                                                                              |
| Calculated MW         | 55,75kDa                                                                                    |
| Storage               | Store at 4°C in dark for 6 months                                                           |

## Application Details

### Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.