

CCL4/MIP1 beta Conjugated Antibody

Catalog No: #C56669

Package Size: #C56669-AF350 100ul #C56669-AF405 100ul #C56669-AF488 100ul #C56669-AF555 100ul #C56669-AF594 100ul
#C56669-AF647 100ul #C56669-AF680 100ul #C56669-AF750 100ul #C56669-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CCL4/MIP1 beta Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse
Specificity	CCL4/MIP1 beta Antibody detects endogenous levels of total CCL4/MIP1 beta
Immunogen Description	A synthesized peptide derived from human CCL4/MIP1 beta
Conjugates	AF Dye
Other Names	ACT2; G-26; HC21; LAG1; LAG-1; MIP1B; SCYA2; SCYA4; MIP1B1; AT744.1;MIP-1-beta;
Accession No.	Uniprot:P13236
Calculated MW	12kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.