

AP2 alpha Conjugated Antibody

Catalog No: #C56689

Package Size: #C56689-AF350 100ul #C56689-AF405 100ul #C56689-AF488 100ul #C56689-AF555 100ul #C56689-AF594 100ul
#C56689-AF647 100ul #C56689-AF680 100ul #C56689-AF750 100ul #C56689-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	AP2 alpha Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	AP2 alpha Antibody detects endogenous levels of total AP2 alpha
Immunogen Description	A synthesized peptide derived from human AP2 alpha
Conjugates	AF Dye
Other Names	AP 2alpha; AP-2; AP2; AP2TF; BOFS; TFAP2; TFAP2A;
Accession No.	Uniprot:P05549
Calculated MW	48kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.