

ADA Conjugated Antibody

Catalog No: #C56699

Package Size: #C56699-AF350 100ul #C56699-AF405 100ul #C56699-AF488 100ul #C56699-AF555 100ul #C56699-AF594 100ul
#C56699-AF647 100ul #C56699-AF680 100ul #C56699-AF750 100ul #C56699-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ADA Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Rat
Specificity	ADA Antibody detects endogenous levels of total ADA
Immunogen Description	A synthesized peptide derived from human ADA
Conjugates	AF Dye
Other Names	ada; ADA1;
Accession No.	Uniprot:P00813
Calculated MW	41kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.