

# Pannexin 1 Conjugated Antibody

Catalog No: #C56718

Package Size: #C56718-AF350 100ul #C56718-AF405 100ul #C56718-AF488 100ul #C56718-AF555 100ul #C56718-AF594 100ul  
#C56718-AF647 100ul #C56718-AF680 100ul #C56718-AF750 100ul #C56718-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                   |
|-----------------------|-------------------------------------------------------------------|
| Product Name          | Pannexin 1 Conjugated Antibody                                    |
| Host Species          | Rabbit                                                            |
| Clonality             | Monoclonal                                                        |
| Isotype               | Rabbit IgG                                                        |
| Purification          | Affinity-chromatography                                           |
| Species Reactivity    | Human                                                             |
| Specificity           | Pannexin 1 Antibody detects endogenous levels of total Pannexin 1 |
| Immunogen Description | A synthesized peptide derived from human Pannexin 1               |
| Conjugates            | AF Dye                                                            |
| Other Names           | innexin; MRS1; Pannexin 1; Panx1; PX1;                            |
| Accession No.         | Uniprot:Q96RD7                                                    |
| Calculated MW         | 48kDa                                                             |
| Storage               | Store at 4°C in dark for 6 months                                 |

## Application Details

### Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.