

Tbx3 Conjugated Antibody

Catalog No: #C56724

Package Size: #C56724-AF350 100ul #C56724-AF405 100ul #C56724-AF488 100ul #C56724-AF555 100ul #C56724-AF594 100ul
#C56724-AF647 100ul #C56724-AF680 100ul #C56724-AF750 100ul #C56724-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tbx3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Tbx3 Antibody detects endogenous levels of total Tbx3
Immunogen Description	A synthesized peptide derived from human Tbx3
Conjugates	AF Dye
Other Names	T box 3; T-box protein 3; Tbx3; TBX3 ISO; UMS; XHL;
Accession No.	Uniprot:O15119
Calculated MW	85kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.