

## GRIK2 Conjugated Antibody

Catalog No: #C56745

Package Size: #C56745-AF350 100ul #C56745-AF405 100ul #C56745-AF488 100ul #C56745-AF555 100ul #C56745-AF594 100ul  
#C56745-AF647 100ul #C56745-AF680 100ul #C56745-AF750 100ul #C56745-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                         |
|-----------------------|---------------------------------------------------------|
| Product Name          | GRIK2 Conjugated Antibody                               |
| Host Species          | Rabbit                                                  |
| Clonality             | Monoclonal                                              |
| Isotype               | Rabbit IgG                                              |
| Purification          | Affinity-chromatography                                 |
| Species Reactivity    | Human Mouse Rat                                         |
| Specificity           | GRIK2 Antibody detects endogenous levels of total GRIK2 |
| Immunogen Description | A synthesized peptide derived from human GRIK2          |
| Conjugates            | AF Dye                                                  |
| Other Names           | EAA4; GLR6; MRT6; GLUK6; GLUR6; GluK2;                  |
| Accession No.         | Uniprot:Q13002                                          |
| Calculated MW         | 103kDa                                                  |
| Storage               | Store at 4°C in dark for 6 months                       |

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.