

M6PR Conjugated Antibody

Catalog No: #C56759



Package Size: #C56759-AF350 100ul #C56759-AF405 100ul #C56759-AF488 100ul #C56759-AF555 100ul #C56759-AF594 100ul
#C56759-AF647 100ul #C56759-AF680 100ul #C56759-AF750 100ul #C56759-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	M6PR Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	M6PR Antibody detects endogenous levels of total M6PR
Immunogen Description	A synthesized peptide derived from human M6PR
Conjugates	AF Dye
Other Names	CD MPR; M6pr; Man6PR; MPR46; MPRD; SMPR;
Accession No.	Uniprot:P20645
Calculated MW	42kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

- AF350 conjugated: most applications: 1: 50 - 1: 250
- AF405 conjugated: most applications: 1: 50 - 1: 250
- AF488 conjugated: most applications: 1: 50 - 1: 250
- AF555 conjugated: most applications: 1: 50 - 1: 250
- AF594 conjugated: most applications: 1: 50 - 1: 250
- AF647 conjugated: most applications: 1: 50 - 1: 250
- AF680 conjugated: most applications: 1: 50 - 1: 250
- AF750 conjugated: most applications: 1: 50 - 1: 250
- Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.