

# CELF1 Conjugated Antibody

Catalog No: #C56800

Package Size: #C56800-AF350 100ul #C56800-AF405 100ul #C56800-AF488 100ul #C56800-AF555 100ul #C56800-AF594 100ul  
#C56800-AF647 100ul #C56800-AF680 100ul #C56800-AF750 100ul #C56800-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                               |
|-----------------------|---------------------------------------------------------------|
| Product Name          | CELF1 Conjugated Antibody                                     |
| Host Species          | Rabbit                                                        |
| Clonality             | Monoclonal                                                    |
| Isotype               | Rabbit IgG                                                    |
| Purification          | Affinity-chromatography                                       |
| Species Reactivity    | Human Mouse Rat                                               |
| Specificity           | CELF1 Antibody detects endogenous levels of total CELF1       |
| Immunogen Description | A synthesized peptide derived from human CELF1                |
| Conjugates            | AF Dye                                                        |
| Other Names           | BRUNOL2; celf1; CUGBP; CUGBP1; EDEN BP; hNab50; NAB50; NAPOR; |
| Accession No.         | Uniprot:Q92879                                                |
| Calculated MW         | 52kDa                                                         |
| Storage               | Store at 4°C in dark for 6 months                             |

## Application Details

### Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.