

TCF4 Conjugated Antibody

Catalog No: #C56801

Package Size: #C56801-AF350 100ul #C56801-AF405 100ul #C56801-AF488 100ul #C56801-AF555 100ul #C56801-AF594 100ul
#C56801-AF647 100ul #C56801-AF680 100ul #C56801-AF750 100ul #C56801-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TCF4 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	TCF4 Antibody detects endogenous levels of total TCF4
Immunogen Description	A synthesized peptide derived from human TCF4
Conjugates	AF Dye
Other Names	E2-2; ITF2; PTHS; SEF2; ITF-2; SEF-2; TCF-4; SEF2-1; SEF2-1A; SEF2-1B; bHLHb19;
Accession No.	Uniprot:P15884
Calculated MW	71,90kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.