

PI 3 Kinase R4 Conjugated Antibody

Catalog No: #C56805

Package Size: #C56805-AF350 100ul #C56805-AF405 100ul #C56805-AF488 100ul #C56805-AF555 100ul #C56805-AF594 100ul #C56805-AF647 100ul #C56805-AF680 100ul #C56805-AF750 100ul #C56805-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PI 3 Kinase R4 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	PI 3 Kinase R4 Antibody detects endogenous levels of total PI 3 Kinase R4
Immunogen Description	A synthesized peptide derived from human PI 3 Kinase R4
Conjugates	AF Dye
Other Names	P150; Phosphoinositide 3 kinase regulatory subunit 4; PI3 kinase p150 subunit; PIK3R4; VPS15;
Accession No.	Uniprot:Q99570
Calculated MW	150kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.