

## ABI2 Conjugated Antibody

Catalog No: #C56814

Package Size: #C56814-AF350 100ul #C56814-AF405 100ul #C56814-AF488 100ul #C56814-AF555 100ul #C56814-AF594 100ul  
#C56814-AF647 100ul #C56814-AF680 100ul #C56814-AF750 100ul #C56814-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                      |
|-----------------------|------------------------------------------------------------------------------------------------------|
| Product Name          | ABI2 Conjugated Antibody                                                                             |
| Host Species          | Rabbit                                                                                               |
| Clonality             | Monoclonal                                                                                           |
| Isotype               | Rabbit IgG                                                                                           |
| Purification          | Affinity-chromatography                                                                              |
| Species Reactivity    | Human Mouse Rat                                                                                      |
| Specificity           | ABI2 Antibody detects endogenous levels of total ABI2                                                |
| Immunogen Description | A synthesized peptide derived from human ABI2                                                        |
| Conjugates            | AF Dye                                                                                               |
| Other Names           | Abelson interactor 2; ABI2; ABI2B; Abl binding protein 3; AblBP3; ArgBP1; ARGBPIA; ArgBPIB; SSH3BP2; |
| Accession No.         | Uniprot:Q9NYB9                                                                                       |
| Calculated MW         | 75kDa                                                                                                |
| Storage               | Store at 4°C in dark for 6 months                                                                    |

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.