

Eph receptor B3 Conjugated Antibody

Catalog No: #C56827

Package Size: #C56827-AF350 100ul #C56827-AF405 100ul #C56827-AF488 100ul #C56827-AF555 100ul #C56827-AF594 100ul
 #C56827-AF647 100ul #C56827-AF680 100ul #C56827-AF750 100ul #C56827-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	Eph receptor B3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Eph receptor B3 Antibody detects endogenous levels of total Eph receptor B3
Immunogen Description	A synthesized peptide derived from human Eph receptor B3
Conjugates	AF Dye
Other Names	Cek10; EK2; Embryonic kinase 2; ephb3; ETK2; hEK2; Mdk5; Sek4; TYRO6;
Accession No.	Uniprot:P54753
Calculated MW	110kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.