

## Lgi1 Conjugated Antibody

Catalog No: #C56831

Package Size: #C56831-AF350 100ul #C56831-AF405 100ul #C56831-AF488 100ul #C56831-AF555 100ul #C56831-AF594 100ul  
#C56831-AF647 100ul #C56831-AF680 100ul #C56831-AF750 100ul #C56831-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                 |
|-----------------------|-----------------------------------------------------------------|
| Product Name          | Lgi1 Conjugated Antibody                                        |
| Host Species          | Rabbit                                                          |
| Clonality             | Monoclonal                                                      |
| Isotype               | Rabbit IgG                                                      |
| Purification          | Affinity-chromatography                                         |
| Species Reactivity    | Human Mouse Rat                                                 |
| Specificity           | Lgi1 Antibody detects endogenous levels of total Lgi1           |
| Immunogen Description | A synthesized peptide derived from human Lgi1                   |
| Conjugates            | AF Dye                                                          |
| Other Names           | ADLTE; ADPAEF; ADPEAF; Epitempin 1; EPITEMPIN; EPT; ETL1; LGI1; |
| Accession No.         | Uniprot:O95970                                                  |
| Calculated MW         | 64kDa                                                           |
| Storage               | Store at 4°C in dark for 6 months                               |

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.