

PEX19 Conjugated Antibody

Catalog No: #C56838

Package Size: #C56838-AF350 100ul #C56838-AF405 100ul #C56838-AF488 100ul #C56838-AF555 100ul #C56838-AF594 100ul
#C56838-AF647 100ul #C56838-AF680 100ul #C56838-AF750 100ul #C56838-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PEX19 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Rat
Specificity	PEX19 Antibody detects endogenous levels of total PEX19
Immunogen Description	A synthesized peptide derived from human PEX19
Conjugates	AF Dye
Other Names	HK33; PBD12A; Peroxin 19; PEX19; PMP1; PMPI; PXF; PXMP1;
Accession No.	Uniprot:P40855
Calculated MW	35,40kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.