

MYH7B Conjugated Antibody

Catalog No: #C56846



Package Size: #C56846-AF350 100ul #C56846-AF405 100ul #C56846-AF488 100ul #C56846-AF555 100ul #C56846-AF594 100ul #C56846-AF647 100ul #C56846-AF680 100ul #C56846-AF750 100ul #C56846-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	MYH7B Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	MYH7B Antibody detects endogenous levels of total MYH7B
Immunogen Description	A synthesized peptide derived from human MYH7B
Conjugates	AF Dye
Other Names	Myh7b; Myosin-7B;
Accession No.	Uniprot:A7E2Y1
Calculated MW	221kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:
AF350 conjugated: most applications: 1: 50 - 1: 250
AF405 conjugated: most applications: 1: 50 - 1: 250
AF488 conjugated: most applications: 1: 50 - 1: 250
AF555 conjugated: most applications: 1: 50 - 1: 250
AF594 conjugated: most applications: 1: 50 - 1: 250
AF647 conjugated: most applications: 1: 50 - 1: 250
AF680 conjugated: most applications: 1: 50 - 1: 250
AF750 conjugated: most applications: 1: 50 - 1: 250
Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.