

CCL23 Conjugated Antibody

Catalog No: #C56887



Package Size: #C56887-AF350 100ul #C56887-AF405 100ul #C56887-AF488 100ul #C56887-AF555 100ul #C56887-AF594 100ul
#C56887-AF647 100ul #C56887-AF680 100ul #C56887-AF750 100ul #C56887-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CCL23 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	CCL23 Antibody detects endogenous levels of total CCL23
Immunogen Description	A synthesized peptide derived from human CCL23
Conjugates	AF Dye
Other Names	CCL23; CKB8; Hmrp-2a; MIP3; MPIF1; SCYA23;
Accession No.	Uniprot:P55773
Calculated MW	13kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB:1:500~1:2000

ICC/IF:1:50~1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.