

PDE4B Conjugated Antibody

Catalog No: #C56888

Package Size: #C56888-AF350 100ul #C56888-AF405 100ul #C56888-AF488 100ul #C56888-AF555 100ul #C56888-AF594 100ul
#C56888-AF647 100ul #C56888-AF680 100ul #C56888-AF750 100ul #C56888-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PDE4B Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse
Specificity	PDE4B Antibody detects endogenous levels of total PDE4B
Immunogen Description	A synthesized peptide derived from human PDE4B
Conjugates	AF Dye
Other Names	DPDE4; PDE 4B; PDE IVB; PDE32; Pde4b; PDE4B5; PDEIVB; Phosphodiesterase 4B;
Accession No.	Uniprot:Q07343
Calculated MW	95kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.