

p66 alpha Conjugated Antibody

Catalog No: #C56903

Package Size: #C56903-AF350 100ul #C56903-AF405 100ul #C56903-AF488 100ul #C56903-AF555 100ul #C56903-AF594 100ul
#C56903-AF647 100ul #C56903-AF680 100ul #C56903-AF750 100ul #C56903-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	p66 alpha Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse
Specificity	p66 alpha Antibody detects endogenous levels of total p66 alpha
Immunogen Description	A synthesized peptide derived from human p66 alpha
Conjugates	AF Dye
Other Names	GATAD2A; Hp66alpha; P66 alpha; p66alpha;
Accession No.	Uniprot:Q86YP4
Calculated MW	65,68kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.