

## NAT10 Conjugated Antibody

Catalog No: #C56904

Package Size: #C56904-AF350 100ul #C56904-AF405 100ul #C56904-AF488 100ul #C56904-AF555 100ul #C56904-AF594 100ul  
#C56904-AF647 100ul #C56904-AF680 100ul #C56904-AF750 100ul #C56904-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                            |
|-----------------------|------------------------------------------------------------|
| Product Name          | NAT10 Conjugated Antibody                                  |
| Host Species          | Rabbit                                                     |
| Clonality             | Monoclonal                                                 |
| Isotype               | Rabbit IgG                                                 |
| Purification          | Affinity-chromatography                                    |
| Species Reactivity    | Human Mouse Rat                                            |
| Specificity           | NAT10 Antibody detects endogenous levels of total NAT10    |
| Immunogen Description | A synthesized peptide derived from human NAT10             |
| Conjugates            | AF Dye                                                     |
| Other Names           | ALP; hALP; N acetyltransferase like protein; NAT10; NET43; |
| Accession No.         | Uniprot:Q9H0A0                                             |
| Calculated MW         | 116kDa                                                     |
| Storage               | Store at 4°C in dark for 6 months                          |

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.