

EHMT2 Conjugated Antibody

Catalog No: #C56905

Package Size: #C56905-AF350 100ul #C56905-AF405 100ul #C56905-AF488 100ul #C56905-AF555 100ul #C56905-AF594 100ul #C56905-AF647 100ul #C56905-AF680 100ul #C56905-AF750 100ul #C56905-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	EHMT2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	EHMT2 Antibody detects endogenous levels of total EHMT2
Immunogen Description	A synthesized peptide derived from EHMT2
Conjugates	AF Dye
Other Names	Bat8; Ehmt2; Euchromatic histone lysine methyltransferase 2; G9a protein; GAT8; Histone H3 K9 methyltransferase 3; NG36;
Accession No.	Uniprot:Q96KQ7
Calculated MW	160,170kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.