

MRRF Conjugated Antibody

Catalog No: #C56965

Package Size: #C56965-AF350 100ul #C56965-AF405 100ul #C56965-AF488 100ul #C56965-AF555 100ul #C56965-AF594 100ul
#C56965-AF647 100ul #C56965-AF680 100ul #C56965-AF750 100ul #C56965-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MRRF Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse
Specificity	MRRF Antibody detects endogenous levels of total MRRF
Immunogen Description	A synthesized peptide derived from human MRRF
Conjugates	AF Dye
Other Names	MRFF; MTRRF; Ribosome recycling factor, mitochondrial; RRF;
Accession No.	Uniprot:Q96E11
Calculated MW	24kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.