

LTB4R Conjugated Antibody

Catalog No: #C56973

Package Size: #C56973-AF350 100ul #C56973-AF405 100ul #C56973-AF488 100ul #C56973-AF555 100ul #C56973-AF594 100ul
#C56973-AF647 100ul #C56973-AF680 100ul #C56973-AF750 100ul #C56973-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LTB4R Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	LTB4R Antibody detects endogenous levels of total LTB4R
Immunogen Description	A synthesized peptide derived from human LTB4R
Conjugates	AF Dye
Other Names	BLT; BLT1; BLTR; CMKRL1; GPR16; LTB4R1; LTBR1; P2RY7; P2Y7;
Accession No.	Uniprot:Q15722
Calculated MW	38kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.