

CD84 Conjugated Antibody

Catalog No: #C56987



Package Size: #C56987-AF350 100ul #C56987-AF405 100ul #C56987-AF488 100ul #C56987-AF555 100ul #C56987-AF594 100ul
#C56987-AF647 100ul #C56987-AF680 100ul #C56987-AF750 100ul #C56987-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CD84 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	CD84 Antibody detects endogenous levels of total CD84
Immunogen Description	A synthesized peptide derived from human CD84
Conjugates	AF Dye
Other Names	CD84; CDw84; hCD84; LY9B; mCD84; SLAMF5;
Accession No.	Uniprot:Q9UIB8
Calculated MW	60-85kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

- AF350 conjugated: most applications: 1: 50 - 1: 250
- AF405 conjugated: most applications: 1: 50 - 1: 250
- AF488 conjugated: most applications: 1: 50 - 1: 250
- AF555 conjugated: most applications: 1: 50 - 1: 250
- AF594 conjugated: most applications: 1: 50 - 1: 250
- AF647 conjugated: most applications: 1: 50 - 1: 250
- AF680 conjugated: most applications: 1: 50 - 1: 250
- AF750 conjugated: most applications: 1: 50 - 1: 250
- Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.