

IL15RA Conjugated Antibody

Catalog No: #C56990

Package Size: #C56990-AF350 100ul #C56990-AF405 100ul #C56990-AF488 100ul #C56990-AF555 100ul #C56990-AF594 100ul
#C56990-AF647 100ul #C56990-AF680 100ul #C56990-AF750 100ul #C56990-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	IL15RA Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	IL15RA Antibody detects endogenous levels of total IL15RA
Immunogen Description	A synthesized peptide derived from human IL15RA
Conjugates	AF Dye
Other Names	CD215; IL 15R alpha; IL-15 receptor subunit alpha; IL-15R-alpha; IL-15RA; Il15ra; Interleukin 15 receptor alpha; sIL-15 receptor subunit alpha; sIL-15R-alpha; sIL-15RA ; Soluble interleukin 15 receptor subunit alpha;
Accession No.	Uniprot:Q13261
Calculated MW	38kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB:1:500~1:2000

ICC/IF:1:50~1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.