

CAP1 Conjugated Antibody

Catalog No: #C56995



Package Size: #C56995-AF350 100ul #C56995-AF405 100ul #C56995-AF488 100ul #C56995-AF555 100ul #C56995-AF594 100ul #C56995-AF647 100ul #C56995-AF680 100ul #C56995-AF750 100ul #C56995-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CAP1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	CAP1 Antibody detects endogenous levels of total CAP1
Immunogen Description	A synthesized peptide derived from human CAP1
Conjugates	AF Dye
Other Names	CAP1; CAP1 PEN;
Accession No.	Uniprot:Q01518
Calculated MW	55kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

- AF350 conjugated: most applications: 1: 50 - 1: 250
- AF405 conjugated: most applications: 1: 50 - 1: 250
- AF488 conjugated: most applications: 1: 50 - 1: 250
- AF555 conjugated: most applications: 1: 50 - 1: 250
- AF594 conjugated: most applications: 1: 50 - 1: 250
- AF647 conjugated: most applications: 1: 50 - 1: 250
- AF680 conjugated: most applications: 1: 50 - 1: 250
- AF750 conjugated: most applications: 1: 50 - 1: 250
- Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.