

MC5 Receptor Conjugated Antibody

Catalog No: #C56996

Package Size: #C56996-AF350 100ul #C56996-AF405 100ul #C56996-AF488 100ul #C56996-AF555 100ul #C56996-AF594 100ul
#C56996-AF647 100ul #C56996-AF680 100ul #C56996-AF750 100ul #C56996-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MC5 Receptor Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	MC5 Receptor Antibody detects endogenous levels of total MC5 Receptor
Immunogen Description	A synthesized peptide derived from human MC5 Receptor
Conjugates	AF Dye
Other Names	MC2; MC5R; Melanocortin 5 receptor;
Accession No.	Uniprot:P33032
Calculated MW	45kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.