

SAMM50 Conjugated Antibody

Catalog No: #C57003

Package Size: #C57003-AF350 100ul #C57003-AF405 100ul #C57003-AF488 100ul #C57003-AF555 100ul #C57003-AF594 100ul
#C57003-AF647 100ul #C57003-AF680 100ul #C57003-AF750 100ul #C57003-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SAMM50 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	SAMM50 Antibody detects endogenous levels of total SAMM50
Immunogen Description	A synthesized peptide derived from human SAMM50
Conjugates	AF Dye
Other Names	CGI 51; OMP85; SAM50; Samm50; TOB55;
Accession No.	Uniprot:Q9Y512
Calculated MW	52kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.