

NPY5R Conjugated Antibody

Catalog No: #C57005



Package Size: #C57005-AF350 100ul #C57005-AF405 100ul #C57005-AF488 100ul #C57005-AF555 100ul #C57005-AF594 100ul
#C57005-AF647 100ul #C57005-AF680 100ul #C57005-AF750 100ul #C57005-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	NPY5R Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	NPY5R Antibody detects endogenous levels of total NPY5R
Immunogen Description	A synthesized peptide derived from human NPY5R
Conjugates	AF Dye
Other Names	NPY 5R; NPY R5; NPY Y5; NPY Y5 receptor; NPY5R; NPYR5; NPY5Y5;Y5 receptor;
Accession No.	Uniprot:Q15761
Calculated MW	51kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

- AF350 conjugated: most applications: 1: 50 - 1: 250
- AF405 conjugated: most applications: 1: 50 - 1: 250
- AF488 conjugated: most applications: 1: 50 - 1: 250
- AF555 conjugated: most applications: 1: 50 - 1: 250
- AF594 conjugated: most applications: 1: 50 - 1: 250
- AF647 conjugated: most applications: 1: 50 - 1: 250
- AF680 conjugated: most applications: 1: 50 - 1: 250
- AF750 conjugated: most applications: 1: 50 - 1: 250
- Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.