

PTPN1/2 Conjugated Antibody

Catalog No: #C57007

Package Size: #C57007-AF350 100ul #C57007-AF405 100ul #C57007-AF488 100ul #C57007-AF555 100ul #C57007-AF594 100ul
#C57007-AF647 100ul #C57007-AF680 100ul #C57007-AF750 100ul #C57007-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PTPN1/2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse
Specificity	PTPN1/2 Antibody detects endogenous levels of total PTPN1/2
Immunogen Description	A synthesized peptide derived from human PTPN1/2
Conjugates	AF Dye
Other Names	Protein tyrosine phosphatase non receptor type 2; PTN2; Ptpn2; PTPT; T cell protein tyrosine phosphatase; TC PTP; TCELLPTP; TCPTP; Tyrosine protein phosphatase non receptor type 2;
Accession No.	Uniprot:P18031/P17706
Calculated MW	40-48kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.